

## ACCLINO<sup>®</sup> flex plus Stent

- » Highly flexible closed cell design
- » More treatment options
- » Compatible with 0.0165" – 0.027" ID microcatheters



# ACCLINO<sup>®</sup> flex plus

## Stent

### Highly flexible

Self-expanding laser-cut Nitinol stent for the treatment of intracranial aneurysms.

### More treatment options

Increased range of sizes, suitable for vessel diameters from 1.5 – 7.0 mm. For an easy handling all sizes up to 5.5 mm in diameter are deliverable through 0.0165" – 0.017" ID microcatheters allowing sequential stent and coil placement without changing the microcatheter.



**SAFE & EFFICIENT<sup>1</sup>**

**100 %**

Technical success

**92.3 %**

Complete occlusion at early follow-up (6 months)


**11 %**

Retreatment

## Flexible

Patented asymmetric closed cell design of the ACCLINO® flex plus Stent leads to highest flexibility in its class of laser-cut stents.

The stent features excellent vessel wall apposition and conformability even in challenging anatomies.




Highest flexibility

## Smooth

Extra flat stent structure and perfectly electropolished surface enable smooth delivery.

Stent structure ensures high permeability in the vessel lumen and leads to a low thrombogenicity.<sup>2</sup>

The ACCLINO® flex plus Stent guarantees reliable coil retention with minimal metal coverage.



Low and smooth profile

## Secure

Easy and clear radiopaque marker concept and good visibility ensure accurate positioning as well as optimal control during procedure.

Repositionable up to 90 % of its length.

Transport wire with Acandis' s.e.c.u.r.e. GP Technology to enhance the grip and push for a controlled and safe stent placement.



Good visibility


## Radiopaque Marker Concept

### Distal transport wire marker

for constant control of the transport wire tip during procedure

### Three distal / proximal device markers

for additional control of stent position and opening behaviour



### Central transport wire marker

for perfect alignment under the aneurysm neck

### Proximal transport wire marker

indicating the „point of no return“ up to which the stent can be repositioned securely

### Distal micro-catheter marker

# ORDERING INFORMATION | ACCLINO® flex plus

Labelled ACCLINO® flex plus Dimensions (mm)	Reference Number	Stent Diameter (mm)	Stent Length (mm)	Recommended Vessel Diameter (mm)	Recommended (Required) Catheters for Delivery (Inch)	
HRF 3.0 × 15 3.0 × 20 3.0 × 25 3.0 × 30 3.0 × 35	01-001122	3.0	15	1.5–2.5	NeuroSlider® 17 NeuroSlider® 17 DLC (ID: 0.0165 – 0.017)	
	01-001123	3.0	20	1.5–2.5		
	01-001124	3.0	25	1.5–2.5		
	01-001125	3.0	30	1.5–2.5		
	01-001126	3.0	35	1.5–2.5		
HRF 3.5 × 15 3.5 × 20 3.5 × 25 3.5 × 30 3.5 × 35	01-001132	3.5	15	1.5–3.0		
	01-001133	3.5	20	1.5–3.0		
	01-001134	3.5	25	1.5–3.0		
	01-001135	3.5	30	1.5–3.0		
	01-001136	3.5	35	1.5–3.0		
HRF 4.0 × 15 4.0 × 20 4.0 × 25 4.0 × 30 4.0 × 35	01-001142	4.0	15	2.5–3.5		
	01-001143	4.0	20	2.5–3.5		
	01-001144	4.0	25	2.5–3.5		
	01-001145	4.0	30	2.5–3.5		
	01-001146	4.0	35	2.5–3.5		
HRF 4.5 × 15 4.5 × 20 4.5 × 25 4.5 × 30 4.5 × 35	01-001152	4.5	15	2.5–4.0		
	01-001153	4.5	20	2.5–4.0		
	01-001154	4.5	25	2.5–4.0		
	01-001155	4.5	30	2.5–4.0		
	01-001156	4.5	35	2.5–4.0		
HRF 5.0 × 15 5.0 × 20 5.0 × 25 5.0 × 30 5.0 × 35	01-001162	5.0	15	3.0–4.5		
	01-001163	5.0	20	3.0–4.5		
	01-001164	5.0	25	3.0–4.5		
	01-001165	5.0	30	3.0–4.5		
	01-001166	5.0	35	3.0–4.5		
HRF 5.5 × 20 5.5 × 25 5.5 × 30 5.5 × 35	01-001173	5.5	20	3.5–5.0		
	01-001174	5.5	25	3.5–5.0		
	01-001175	5.5	30	3.5–5.0		
	01-001176	5.5	35	3.5–5.0		
	HRF 6.5 × 20 6.5 × 25 6.5 × 30 6.5 × 35	01-001193	6.5	20		4.0–6.0
01-001194		6.5	25	4.0–6.0		
01-001195		6.5	30	4.0–6.0		
01-001196		6.5	35	4.0–6.0		
HRF 8.0 × 20 8.0 × 30 8.0 × 40 8.0 × 60	01-001213	8.0	20	6.0–7.0		NeuroSlider® 27 NeuroSlider® 27 DLC (ID: 0.027)
	01-001215	8.0	30	6.0–7.0		
	01-001217	8.0	40	6.0–7.0		
	01-001221	8.0	60	6.0–7.0		

**HRF:** High Radial Force – compared to ACCLINO® flex plus Stents within the same recommended vessel diameter.

All changes or modifications, may they be technical or other, or changes in the availability of products are expressly reserved.

Please contact your local Acandis® representative for product availability and information on compatible (micro)catheters.

Not available for sale in the United States.

CE 0297