


## **NeuroSlider<sup>®</sup> DLC** (Micro)Catheter


- » Superior torqueability and pushability
- » Smooth and safe device delivery
- » Long-lasting tip shape retention




# ORDERING INFORMATION | NeuroSlider® DLC

Product Name	Reference Number	ID (Inch)	OD dist. / prox. (French)	Usable Length (cm)	Tip Shape	Tip Marker
NeuroSlider® 17 DLC	01-000282	0.0165	1.9 / 2.3	155	Straight (shapeable)	2
	01-000283	0.0165	1.9 / 2.3	160	Straight (shapeable)	2
	01-000284	0.0165	1.9 / 2.3	167	Straight (shapeable)	2
NeuroSlider® 21 DLC	01-000292	0.021	2.3 / 2.6	155	Straight (shapeable)	2
	01-000293	0.021	2.3 / 2.6	160	Straight (shapeable)	2
	01-000294	0.021	2.3 / 2.6	167	Straight (shapeable)	2
NeuroSlider® 27 DLC	01-000276	0.027	3.0 / 3.1	155	Straight (shapeable)	1
NeuroSlider® 39 DLC	01-000262	0.039	4.0 / 4.1	125	Straight (shapeable)	1
	01-000263	0.039	4.0 / 4.1	135	Straight (shapeable)	1
	01-000264	0.039	4.0 / 4.1	145	Straight (shapeable)	1
NeuroSlider® 52 DLC	01-000252	0.052	5.0 / 5.1	105	Straight (shapeable)	1
	01-000253	0.052	5.0 / 5.1	115	Straight (shapeable)	1
	01-000254	0.052	5.0 / 5.1	125	Straight (shapeable)	1
	01-000255	0.052	5.0 / 5.1	135	Straight (shapeable)	1
	01-000256	0.052	5.0 / 5.1	145	Straight (shapeable)	1

## NeuroSlider® DLC (with 2 tip markers)



## NeuroSlider® DLC (with 1 tip marker)



All changes or modifications, may they be technical or other, or changes in the availability of products are expressly reserved.  
 Please contact your local Acandis® representative for product availability and information on compatible devices.  
 Not available for sale in the United States.