

# EXTesia

## Introducer Sheath Set With Spring Guide Wire/Metallic Needle

Each set contains a sheath, a dilator, a spring guide wire, a metallic needle with gauge compatible with guide wire size selected, a 2.5ml syringe and a scalpel



WWW.CURATIAMED.COM

# EXTesia Introducer Sheath Set

Available in 4,5,6 and 7F sizes with 5,7,10,16 and 25cm in length  
Designed to reduce radial access trauma and spasm

- Sheath material with kink recover technology for maintain its atraumatic shape during procedure & reduce the need to exchange the sheath
- Smooth transition, taper tip and hydrophilic coating to reduce insertion force
- Variety selection of guide wire sizes in 0.018, 0.021, 0.025, 0.035 and 0.038"
- 2 types of needle: metallic and IV catheter gage matching variety of guide wire sizes
- Hemostasis valve with back-bleeding prevention features



Sheath			Guide Wire						
F Size	Working Length (cm)	Length (cm)	Tip Shape	OD (inch)					
				0.018	0.021	0.025	0.035	0.038	
4F	5	45	Straight	2IS040518SSN	2IS040521SSN	2IS040525SSN	2IS040535SSN	-	
			J	2IS040518SJN	2IS040521SJN	2IS040525SJN	2IS040535SJN	-	
	7		Straight	2IS040718SSN	2IS040721SSN	2IS040725SSN	2IS040735SSN	-	
			J	2IS040718SJN	2IS040721SJN	2IS040725SJN	2IS040735SJN	-	
	10		Straight	2IS041018SSN	2IS041021SSN	2IS041025SSN	2IS041035SSN	-	
			J	2IS041018SJN	2IS041021SJN	2IS041025SJN	2IS041035SJN	-	
	16		80	Straight	2IS041618SSN	2IS041621SSN	2IS041625SSN	2IS041635SSN	-
				J	2IS041618SJN	2IS041621SJN	2IS041625SJN	2IS041635SJN	-
	25			Straight	2IS042518SSN	2IS042521SSN	2IS042525SSN	2IS042535SSN	-
				J	2IS042518SJN	2IS042521SJN	2IS042525SJN	2IS042535SJN	-
5F	5	45		Straight	2IS050518SSN	2IS050521SSN	2IS050525SSN	2IS050535SSN	2IS050538SSN
				J	2IS050518SJN	2IS050521SJN	2IS050525SJN	2IS050535SJN	2IS050538SJN
	7			Straight	2IS050718SSN	2IS050721SSN	2IS050725SSN	2IS050735SSN	2IS050738SSN
				J	2IS050718SJN	2IS050721SJN	2IS050725SJN	2IS050735SJN	2IS050738SJN
	10			Straight	2IS051018SSN	2IS051021SSN	2IS051025SSN	2IS051035SSN	2IS051038SSN
				J	2IS051018SJN	2IS051021SJN	2IS051025SJN	2IS051035SJN	2IS051038SJN
	16		80	Straight	2IS051618SSN	2IS051621SSN	2IS051625SSN	2IS051635SSN	2IS051638SSN
				J	2IS051618SJN	2IS051621SJN	2IS051625SJN	2IS051635SJN	2IS051638SJN
	25			Straight	2IS052518SSN	2IS052521SSN	2IS052525SSN	2IS052535SSN	2IS052538SSN
				J	2IS052518SJN	2IS052521SJN	2IS052525SJN	2IS052535SJN	2IS052538SJN
6F	5	45		Straight	2IS060518SSN	2IS060521SSN	2IS060525SSN	2IS060535SSN	2IS060538SSN
				J	2IS060518SJN	2IS060521SJN	2IS060525SJN	2IS060535SJN	2IS060538SJN
	7			Straight	2IS060718SSN	2IS060721SSN	2IS060725SSN	2IS060735SSN	2IS060738SSN
				J	2IS060718SJN	2IS060721SJN	2IS060725SJN	2IS060735SJN	2IS060738SJN
	10			Straight	2IS061018SSN	2IS061021SSN	2IS061025SSN	2IS061035SSN	2IS061038SSN
				J	2IS061018SJN	2IS061021SJN	2IS061025SJN	2IS061035SJN	2IS061038SJN
	16		80	Straight	2IS061618SSN	2IS061621SSN	2IS061625SSN	2IS061635SSN	2IS061638SSN
				J	2IS061618SJN	2IS061621SJN	2IS061625SJN	2IS061635SJN	2IS061638SJN
	25			Straight	2IS062518SSN	2IS062521SSN	2IS062525SSN	2IS062535SSN	2IS062538SSN
				J	2IS062518SJN	2IS062521SJN	2IS062525SJN	2IS062535SJN	2IS062538SJN
7F	5	45		Straight	2IS070518SSN	2IS070521SSN	2IS070525SSN	2IS070535SSN	2IS070538SSN
				J	2IS070518SJN	2IS070521SJN	2IS070525SJN	2IS070535SJN	2IS070538SJN
	7			Straight	2IS070718SSN	2IS070721SSN	2IS070725SSN	2IS070735SSN	2IS070738SSN
				J	2IS070718SJN	2IS070721SJN	2IS070725SJN	2IS070735SJN	2IS070738SJN
	10			Straight	2IS071018SSN	2IS071021SSN	2IS071025SSN	2IS071035SSN	2IS071038SSN
				J	2IS071018SJN	2IS071021SJN	2IS071025SJN	2IS071035SJN	2IS071038SJN
	16		80	Straight	2IS071618SSN	2IS071621SSN	2IS071625SSN	2IS071635SSN	2IS071638SSN
				J	2IS071618SJN	2IS071621SJN	2IS071625SJN	2IS071635SJN	2IS071638SJN
	25			Straight	2IS072518SSN	2IS072521SSN	2IS072525SSN	2IS072535SSN	2IS072538SSN
				J	2IS072518SJN	2IS072521SJN	2IS072525SJN	2IS072535SJN	2IS072538SJN

CURATIA MEDICAL

Web: [www.curatiamed.com](http://www.curatiamed.com) E-mail: [info@curatiamed.com](mailto:info@curatiamed.com)

CHINA

Add: 198 Xiangjiang Rd, New District, Suzhou, Jiangsu

UNITED STATES

Add: 830 Stewart Drive, Sunnyvale, CA94085

GERMANY


Add: Blasewitzer Str.41,01307 Dresden

# EXTesia

## Introducer Sheath Set

### With Polymer Jacket Guide Wire/ I.V.Catheter

Each set contains a sheath, a dilator, a polymer jacket guide wire, an I.V.Catheter with gage compatible with guide wire size selected, a 2.5ml syringe and a scalpel




Sheath			Guide Wire				
F Size	Working Length (cm)	Length (cm)	Tip Shape	OD (inch)			
				0.021	0.025	0.035	0.038
4F	5	45	Straight	2IS040521PSI	2IS040525PSI	2IS040535PSI	-
			Angled	2IS040521PAI	2IS040525PAI	2IS040535PAI	-
			Straight	2IS040721PSI	2IS040725PSI	2IS040735PSI	-
	Angled		2IS040721PAI	2IS040725PAI	2IS040735PAI	-	
	Straight		2IS041021PSI	2IS041025PSI	2IS041035PSI	-	
	Angled		2IS041021PAI	2IS041025PAI	2IS041035PAI	-	
	10	80	Straight	2IS041621PSI	2IS041625PSI	2IS041635PSI	-
			Angled	2IS041621PAI	2IS041625PAI	2IS041635PAI	-
			Straight	2IS042521PSI	2IS042525PSI	2IS042535PSI	-
		Angled	2IS042521PAI	2IS042525PAI	2IS042535PAI	-	
		Straight	2IS050521PSI	2IS050525PSI	2IS050535PSI	2IS050538PSI	
		Angled	2IS050521PAI	2IS050525PAI	2IS050535PAI	2IS050538PAI	
5F	5	45	Straight	2IS050721PSI	2IS050725PSI	2IS050735PSI	2IS050738PSI
			Angled	2IS050721PAI	2IS050725PAI	2IS050735PAI	2IS050738PAI
			Straight	2IS051021PSI	2IS051025PSI	2IS051035PSI	2IS051038PSI
	Angled		2IS051021PAI	2IS051025PAI	2IS051035PAI	2IS051038PAI	
	Straight		2IS051621PSI	2IS051625PSI	2IS051635PSI	2IS051638PSI	
	Angled		2IS051621PAI	2IS051625PAI	2IS051635PAI	2IS051638PAI	
	10	80	Straight	2IS052521PSI	2IS052525PSI	2IS052535PSI	2IS052538PSI
			Angled	2IS052521PAI	2IS052525PAI	2IS052535PAI	2IS052538PAI
			Straight	2IS060521PSI	2IS060525PSI	2IS060535PSI	2IS060538PSI
		Angled	2IS060521PAI	2IS060525PAI	2IS060535PAI	2IS060538PAI	
		Straight	2IS060721PSI	2IS060725PSI	2IS060735PSI	2IS060738PSI	
		Angled	2IS060721PAI	2IS060725PAI	2IS060735PAI	2IS060738PAI	
6F	5	45	Straight	2IS061021PSI	2IS061025PSI	2IS061035PSI	2IS061038PSI
			Angled	2IS061021PAI	2IS061025PAI	2IS061035PAI	2IS061038PAI
			Straight	2IS061621PSI	2IS061625PSI	2IS061635PSI	2IS061638PSI
	Angled		2IS061621PAI	2IS061625PAI	2IS061635PAI	2IS061638PAI	
	Straight		2IS062521PSI	2IS062525PSI	2IS062535PSI	2IS062538PSI	
	Angled		2IS062521PAI	2IS062525PAI	2IS062535PAI	2IS062538PAI	
	10	80	Straight	2IS070521PSI	2IS070525PSI	2IS070535PSI	2IS070538PSI
			Angled	2IS070521PAI	2IS070525PAI	2IS070535PAI	2IS070538PSI
			Straight	2IS070721PSI	2IS070725PSI	2IS070735PSI	2IS070738PSI
		Angled	2IS070721PAI	2IS070725PAI	2IS070735PAI	2IS070738PSI	
		Straight	2IS071021PSI	2IS071025PSI	2IS071035PSI	2IS071038PSI	
		Angled	2IS071021PAI	2IS071025PAI	2IS071035PAI	2IS071038PSI	
7F	5	45	Straight	2IS071621PSI	2IS071625PSI	2IS071635PSI	2IS071638PSI
			Angled	2IS071621PAI	2IS071625PAI	2IS071635PAI	2IS071638PSI
			Straight	2IS072521PSI	2IS072525PSI	2IS072535PSI	2IS072538PSI
	Angled		2IS072521PAI	2IS072525PAI	2IS072535PAI	2IS072538PSI	
	Straight		2IS070721PSI	2IS070725PSI	2IS070735PSI	2IS070738PSI	
	Angled		2IS070721PAI	2IS070725PAI	2IS070735PAI	2IS070738PSI	
	10	80	Straight	2IS071021PSI	2IS071025PSI	2IS071035PSI	2IS071038PSI
			Angled	2IS071021PAI	2IS071025PAI	2IS071035PAI	2IS071038PSI
			Straight	2IS071621PSI	2IS071625PSI	2IS071635PSI	2IS071638PSI
		Angled	2IS071621PAI	2IS071625PAI	2IS071635PAI	2IS071638PSI	
		Straight	2IS072521PSI	2IS072525PSI	2IS072535PSI	2IS072538PSI	
		Angled	2IS072521PAI	2IS072525PAI	2IS072535PAI	2IS072538PSI	

# EXTesia

## Introducer Sheath Set

### With Spring Guide Wire/ I.V.Catheter

Each set contains a sheath, a dilator, a spring guide wire, an I.V.Catheter w. gauge compatible guide wire size selected, a 2.5ml syringe and a scalpel



Sheath			Guide Wire					
F Size	Working Length (cm)	Length (cm)	Tip Shape	OD (inch)				
				0.018	0.021	0.025	0.035	0.038
4F	5	45	Straight	2IS040518SSI	2IS040521SSI	2IS040525SSI	2IS040535SSI	-
			J	2IS040518SJI	2IS040521SJI	2IS040525SJI	2IS040535SJI	-
			Straight	2IS040718SSI	2IS040721SSI	2IS040725SSI	2IS040735SSI	-
	J		2IS040718SJI	2IS040721SJI	2IS040725SJI	2IS040735SJI	-	
	Straight		2IS041018SSI	2IS041021SSI	2IS041025SSI	2IS041035SSI	-	
	J		2IS041018SJI	2IS041021SJI	2IS041025SJI	2IS041035SJI	-	
	10	80	Straight	2IS041618SSI	2IS041621SSI	2IS041625SSI	2IS041635SSI	-
			J	2IS041618SJI	2IS041621SJI	2IS041625SJI	2IS041635SJI	-
			Straight	2IS042518SSI	2IS042521SSI	2IS042525SSI	2IS042535SSI	-
		J	2IS042518SJI	2IS042521SJI	2IS042525SJI	2IS042535SJI	-	
		Straight	2IS050518SSI	2IS050521SSI	2IS050525SSI	2IS050535SSI	2IS050538SSI	
		J	2IS050518SJI	2IS050521SJI	2IS050525SJI	2IS050535SJI	2IS050538SJI	
5F	5	45	Straight	2IS050718SSI	2IS050721SSI	2IS050725SSI	2IS050735SSI	2IS050738SSI
			J	2IS050718SJI	2IS050721SJI	2IS050725SJI	2IS050735SJI	2IS050738SJI
			Straight	2IS051018SSI	2IS051021SSI	2IS051025SSI	2IS051035SSI	2IS051038SSI
	J		2IS051018SJI	2IS051021SJI	2IS051025SJI	2IS051035SJI	2IS051038SJI	
	Straight		2IS051618SSI	2IS051621SSI	2IS051625SSI	2IS051635SSI	2IS051638SSI	
	J		2IS051618SJI	2IS051621SJI	2IS051625SJI	2IS051635SJI	2IS051638SJI	
	10	80	Straight	2IS052518SSI	2IS052521SSI	2IS052525SSI	2IS052535SSI	2IS052538SSI
			J	2IS052518SJI	2IS052521SJI	2IS052525SJI	2IS052535SJI	2IS052538SJI
			Straight	2IS060518SSI	2IS060521SSI	2IS060525SSI	2IS060535SSI	2IS060538SSI
		J	2IS060518SJI	2IS060521SJI	2IS060525SJI	2IS060535SJI	2IS060538SJI	
		Straight	2IS060718SSI	2IS060721SSI	2IS060725SSI	2IS060735SSI	2IS060738SSI	
		J	2IS060718SJI	2IS060721SJI	2IS060725SJI	2IS060735SJI	2IS060738SJI	
6F	5	45	Straight	2IS061018SSI	2IS061021SSI	2IS061025SSI	2IS061035SSI	2IS061038SSI
			J	2IS061018SJI	2IS061021SJI	2IS061025SJI	2IS061035SJI	2IS061038SJI
			Straight	2IS061618SSI	2IS061621SSI	2IS061625SSI	2IS061635SSI	2IS061638SSI
	J		2IS061618SJI	2IS061621SJI	2IS061625SJI	2IS061635SJI	2IS061638SJI	
	Straight		2IS062518SSI	2IS062521SSI	2IS062525SSI	2IS062535SSI	2IS062538SSI	
	J		2IS062518SJI	2IS062521SJI	2IS062525SJI	2IS062535SJI	2IS062538SJI	
	10	80	Straight	2IS070518SSI	2IS070521SSI	2IS070525SSI	2IS070535SSI	2IS070538SSI
			J	2IS070518SJI	2IS070521SJI	2IS070525SJI	2IS070535SJI	2IS070538SJI
			Straight	2IS070718SSI	2IS070721SSI	2IS070725SSI	2IS070735SSI	2IS070738SSI
		J	2IS070718SJI	2IS070721SJI	2IS070725SJI	2IS070735SJI	2IS070738SJI	
		Straight	2IS071018SSI	2IS071021SSI	2IS071025SSI	2IS071035SSI	2IS071038SSI	
		J	2IS071018SJI	2IS071021SJI	2IS071025SJI	2IS071035SJI	2IS071038SJI	
7F	5	45	Straight	2IS071618SSI	2IS071621SSI	2IS071625SSI	2IS071635SSI	2IS071638SSI
			J	2IS071618SJI	2IS071621SJI	2IS071625SJI	2IS071635SJI	2IS071638SJI
			Straight	2IS072518SSI	2IS072521SSI	2IS072525SSI	2IS072535SSI	2IS072538SSI
	J		2IS072518SJI	2IS072521SJI	2IS072525SJI	2IS072535SJI	2IS072538SJI	
	Straight		2IS070718SSI	2IS070721SSI	2IS070725SSI	2IS070735SSI	2IS070738SSI	
	J		2IS070718SJI	2IS070721SJI	2IS070725SJI	2IS070735SJI	2IS070738SJI	
	10	80	Straight	2IS071018SSI	2IS071021SSI	2IS071025SSI	2IS071035SSI	2IS071038SSI
			J	2IS071018SJI	2IS071021SJI	2IS071025SJI	2IS071035SJI	2IS071038SJI
			Straight	2IS071618SSI	2IS071621SSI	2IS071625SSI	2IS071635SSI	2IS071638SSI
		J	2IS071618SJI	2IS071621SJI	2IS071625SJI	2IS071635SJI	2IS071638SJI	
		Straight	2IS072518SSI	2IS072521SSI	2IS072525SSI	2IS072535SSI	2IS072538SSI	
		J	2IS072518SJI	2IS072521SJI	2IS072525SJI	2IS072535SJI	2IS072538SJI	