



WWW.CURATIAMED.COM



**XCESS**<sup>®</sup>  
Guiding Catheter

# XCESS<sup>®</sup>

Guiding Catheter

Providing a controlled access



The **XCESS** guiding catheter is designed to provide a pathway through which therapeutic devices are introduced in coronary and peripheral vascular system



- Hybrid braid design yields efficient 1:1 torque control and kink resistance
- Unique large lumen shaft platform enables back-up support while maintaining optimal lumen integrity
- Curve transitions provides flexibility in seating with gentle soft tip engagement
- PTFE liner reduces friction and enhances device delivery



WWW.CURATIAMED.COM

No.	Catheter Tip Shape	Catheter Size(F)	
		5	6
1	Amplatz left 1	102A5L1	102A6L1
2	Amplatz left 2	102A5L2	102A6L2
3	Amplatz left 3	102A5L3	102A6L3
4	Amplatz Right 1	102A5R1	102A6R1
5	Amplatz Right 2	102A5R2	102A6R2
6	Amplatz Right 3	102A5R3	102A6R3
7	Judkins Left 3.0	10215L2	10216L2
8	Judkins Left 3.5	10215L3	10216L3
9	Judkins Left 4.0	10215L4	10216L4
10	Judkins Left 4.5	10215L5	10216L5
11	Judkins Left 5.0	10215L6	10216L6
12	Judkins Left 6.0	10215L7	10216L7
13	Judkins Right 3.0	10215R2	10216R2
14	Judkins Right 3.5	10215R3	10216R3
15	Judkins Right 4.0	10215R4	10216R4
16	Judkins Right 5.0	10215R6	10216R6
17	Judkins Right 6.0	10215R7	10216R7
18	XB 3	102X502	102X602
19	XB 3.5	102X503	102X603
20	XB 4	102X504	102X604
21	XB 4.5	102X505	102X605
22	XBLAD 3	102X5D2	102X6D2
23	XBLAD 3.5	102X5D3	102X6D3

No.	Catheter Tip Shape	Catheter Size(F)	
		5	6
24	XBLAD 4	102X5D4	102X6D4
25	XBLAD 4.5	102X5D5	102X6D5
26	XBC 3	102X5C2	102X6C2
27	XBC 3.5	102X5C3	102X6C3
28	XBC 4	102X5C4	102X6C4
29	XBC 4.5	102X5C5	102X6C5
30	XBRCA	102X5CA	102X6CA
31	XBR 1	102X5R1	102X6R1
32	XBR 2	102X5R2	102X6R2
33	Judkins Curved Left 3.5	102J5L3	102J6L3
34	Judkins Curved Left 4.0	102J5L4	102J6L4
35	Judkins Curved Right 4.0	102J5R4	102J6R4
36	3DRC	10215RC	10216RC
37	Castillo 2	10215C2	10216C2
38	Barbeau	10215BA	10216BA
39	Coronary Bypass-left	102C5L0	102C6L0
40	Coronary Bypass-right	102C5R0	102C6R0
41	El Gamal Bypass	102E500	102E600
42	SON 1	102S501	102S601
43	Internal mammary	102I500	102I600
44	Internal mammary VB-1	102I5VB	102I6VB
45	Multipurpose A	102M5A0	102M6A0
46	Multipurpose B	102M5B0	102M6B0

No.	Catheter Tip Shape	Catheter Size(F)	
		5	6
47	Hockey stick	10215HS	10216HS
48	Modified cerebral burke	10215BL	10216BL
49	Radial brachial	10215RB	10216RB
50	Renal double bumper	102R500	102R600
51	J curve	102150J	102160J

### Catheter Inner Diameter Competitiveness

	Curatia	MDT	Cordis	BSX	Terumo	Merit
	Xcess	Launcher	Vista Brite Tip	Mach1/Convey	Heartrial II	Concierge
5F	0.058	0.058	0.056	0.058	0.058	0.057
6F	0.071	0.071	0.070	0.070/0.071	0.071	0.070

 Curatia Medical Ltd.  
198 Xiangjiang Road, New District, Suzhou, China  
Email: info@curatiamedical.com

Curatia Medical Co.  
830 Stewart Drive, Sunnyvale, CA 94085, USA

Curatia Medical Devices UG  
Blasewitzer Str. 41,01307 Dresden, Germany

